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WORK DETAILS

Engineer Name: SHUYA MIYAZAKI

Event Type: Predictive / Preventative Maintenance

Device Status: Conditional Pass

Service Activity: Preventative Maintenance

EQUIPMENT DETAILS

Product Name: EPIQ 7 Ultrasound System

Serial Number: USO15B0925

Product Number: 795200

Installed Product: 67704928

UDI:

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

Planned Maintenance

System Information

Software Version					7.0.5	
Hardware Version					B.1	

Inspection & Cleaning

Query Customer Performance Issues					Pass	
Inspect Exterior Periph. & Power Cords.					Pass	
Inspect & Clean Trackball					Pass	
Inspect & Clean Monitor					Pass	
Inspect/Clean Control Panels & Keyboard					Pass	
Inspected and cleaned the transducer connector assembly (transducer and system side)?					Yes	
Footswitch, Key- & Remote Control Pad					Not Installed	
Inspect Casters & Brakes Operation					Pass	
Clean/Check System Filters & Fans					Pass	
Check Articulation motion.(monitor, etc)					Fail	ガスストラットの劣化が見られます。交換を推奨いたします。
Check System Battery					Pass	

Performance & Peripherals

Execute Full Set of System Diagnostics					Pass	
Check, Download & Clear Logs					Pass	
Perform Hard-disk Maintenance					Pass	
Backup of System Settings and Configs.					Pass	

Min Max UoM Measurement Result

Comments

Planned Maintenance

Performance & Peripherals

Check Operation All User Ctrls. & Lights					Pass
Check ECG Cable and Trace on Screen					Pass
Check Operation of Major Applications					Pass
Test Printers					Pass
Check Operation Dig. Recording Devices					Pass
Confirm remote connectivity					Not Installed

System Safety General

Protective Earth resistance (<= 200 mOhm)	200	mOhm	30		Pass
Earth Leakage Current. NC, Highest val. (<= 500 µA)	500	µA	193.5		Pass
Earth Leakage Current. SF. Highest val. (<= 1000 µA)	1000	µA	371		Pass

ECG Body - BF

ECG Lead Leakage Current "Norm Cond"	100	µA	0.1		Pass
ECG LL Curr. Highest SF or 2SF reading	500	µA	1.8		Pass
ECG Lead Isolation, ISO/MAP "Sink"	500	µA	1.8		Pass

Transducer Body - BF

Transducer model					L12-3
Transducer serial number					B1LVYF
Transducer Cables, Connectors & Lens					PASS
XDR Leakage Current "Normal Condition" (<= 100 µA)	100		0.1		Pass
XDR Leak Curr. Highest SF or 2SF reading" (<= 500 µA)	500		1.1		Pass
XDR Isolation, ISO/MAP "Sink" (<= 500 µA)	500		3		Pass
Run Transducer Test					PASS

Transducer Body - BF

Transducer model					C9-2
Transducer serial number					B1KXPJ
Transducer Cables, Connectors & Lens					PASS
XDR Leakage Current "Normal Condition" (<= 100 µA)	100		0.1		Pass
XDR Leak Curr. Highest SF or 2SF reading" (<= 500 µA)	500		2.4		Pass
XDR Isolation, ISO/MAP "Sink" (<= 500 µA)	500		5.4		Pass
Run Transducer Test					PASS

Transducer Heart - CF

Transducer model					S9-2
Transducer serial number					B3G9RZ
Transducer Cables, Connectors & Lens					PASS
XDR Leakage Current "Normal Condition" (<= 10 µA)	10		0.1		Pass

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

Planned Maintenance

Transducer Heart - CF

XDR Leak Curr. Highest SF or 2SF reading (<= 50 µA)		50		1.7	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 50 µA)		50		3.7	Pass	
Run Transducer Test					PASS	

Transducer Heart - CF

Transducer model					S8-3	
Transducer serial number					B1JVCG	
Transducer Cables, Connectors & Lens					PASS	
XDR Leakage Current "Normal Condition" (<= 10 µA)		10		0.1	Pass	
XDR Leak Curr. Highest SF or 2SF reading (<= 50 µA)		50		1.7	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 50 µA)		50		1.8	Pass	
Run Transducer Test					PASS	

Transducer Heart - CF

Transducer model					X5-1	
Transducer serial number					F03GNG	
Transducer Cables, Connectors & Lens					PASS	
XDR Leakage Current "Normal Condition" (<= 10 µA)		10		0.1	Pass	
XDR Leak Curr. Highest SF or 2SF reading (<= 50 µA)		50		1.4	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 50 µA)		50		2.9	Pass	
Run Transducer Test					PASS	

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Debrief

Number of Soft Failures 1

Philips Representative Signature:

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