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### WORK DETAILS

Engineer Name: SHUYA MIYAZAKI

Event Type: Predictive / Preventative Maintenance

Device Status: Conditional Pass

Service Activity: Preventative Maintenance

### EQUIPMENT DETAILS

Product Name: EPIQ 7 Ultrasound System

Serial Number: USO15B0926

Product Number: 795200

Installed Product: 67704962

UDI:

	Min	Max	UoM	Measurement	Result	Comments
--	-----	-----	-----	-------------	--------	----------

## Planned Maintenance

### System Information

Software Version					7.0.5	
Hardware Version					B.1	

### Inspection & Cleaning

Query Customer Performance Issues					Pass	
Inspect Exterior Periph. & Power Cords.					Pass	
Inspect & Clean Trackball					Pass	
Inspect & Clean Monitor					Pass	
Inspect/Clean Control Panels & Keyboard					Pass	
Inspected and cleaned the transducer connector assembly (transducer and system side)?					Yes	
Footswitch, Key- & Remote Control Pad					Not Installed	
Inspect Casters & Brakes Operation					Pass	
Clean/Check System Filters & Fans					Pass	
Check Articulation motion.(monitor, etc)					Pass	
Check System Battery					Pass	

### Performance & Peripherals

Execute Full Set of System Diagnostics					Pass	
Check, Download & Clear Logs					Pass	
Perform Hard-disk Maintenance					Pass	
Backup of System Settings and Configs.					Pass	

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

## Planned Maintenance

### Performance & Peripherals

Check Operation All User Ctrls. & Lights				Pass	
Check ECG Cable and Trace on Screen				Pass	
Check Operation of Major Applications				Pass	
Test Printers				Pass	
Check Operation Dig. Recording Devices				Pass	
Confirm remote connectivity				Not Installed	

### System Safety General

Protective Earth resistance (<= 200 mOhm)	200	mOhm	19	Pass	
Earth Leakage Current. NC, Highest val. (<= 500 µA)	500	µA	177.4	Pass	
Earth Leakage Current. SF. Highest val. (<= 1000 µA)	1000	µA	343	Pass	

### ECG Body - BF

ECG Lead Leakage Current "Norm Cond"	100	µA	0.1	Pass	
ECG LL Curr. Highest SF or 2SF reading	500	µA	1.6	Pass	
ECG Lead Isolation, ISO/MAP "Sink"	500	µA	1.8	Pass	

### Transducer Body - BF

Transducer model				L12-3	
Transducer serial number				B1LWB6	
Transducer Cables, Connectors & Lens				PASS	
XDR Leakage Current "Normal Condition" (<= 100 µA)	100		0.1	Pass	
XDR Leak Curr. Highest SF or 2SF reading" (<= 500 µA)	500		1.2	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 500 µA)	500		2.3	Pass	
Run Transducer Test				PASS	

### Transducer Body - BF

Transducer model				C9-2	
Transducer serial number				B1KF0B	
Transducer Cables, Connectors & Lens				PASS	
XDR Leakage Current "Normal Condition" (<= 100 µA)	100		0.1	Pass	
XDR Leak Curr. Highest SF or 2SF reading" (<= 500 µA)	500		2.1	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 500 µA)	500		3.8	Pass	
Run Transducer Test				PASS	

### Transducer Heart - CF

Transducer model				S922	
Transducer serial number				B3LN2R	
Transducer Cables, Connectors & Lens				PASS	
XDR Leakage Current "Normal Condition" (<= 10 µA)	10		0.1	Pass	

	Min	Max	UoM	Measurement	Result	Comments
<b>Planned Maintenance</b>						
<b>Transducer Heart - CF</b>						
XDR Leak Curr. Highest SF or 2SF reading (<= 50 µA)		50		1.6	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 50 µA)		50		2.7	Pass	
Run Transducer Test					PASS	
<b>Transducer Heart - CF</b>						
Transducer model					S8-3	
Transducer serial number					B1JVK5	
Transducer Cables, Connectors & Lens					PASS	
XDR Leakage Current "Normal Condition" (<= 10 µA)		10		0.1	Pass	
XDR Leak Curr. Highest SF or 2SF reading (<= 50 µA)		50		0.8	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 50 µA)		50		1	Pass	
Run Transducer Test					PASS	
<b>Transducer Heart - CF</b>						
Transducer model					X5-1	
Transducer serial number					B3FY20	
Transducer Cables, Connectors & Lens					PASS	
XDR Leakage Current "Normal Condition" (<= 10 µA)		10		0.1	Pass	
XDR Leak Curr. Highest SF or 2SF reading (<= 50 µA)		50		1.6	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 50 µA)		50		2.7	Pass	
Run Transducer Test					PASS	
<b>Transducer TEE - BF</b>						
Transducer model					X8-2t	
Transducer serial number					B38VT7	
Inspect Transducer Cable					PASS	
XDR Leakage Current "Normal Condition" (<= 300 µA)		300	µA	0.1	Pass	
XDR Leakage Current "Single Condition" (<= 500 µA)		500	µA	12.3	Pass	
XDR Isolation, ISO/MAP "Sink" (<= 500 µA)		500	µA	29.3	Pass	
Inspect TEE Xdcrs & Controls					PASS	
Check Sheath condition					PASS	
Check Nose Piece for Damage					PASS	
Check Articulation Deflection Cover					PASS	
Check Handle / Connector for Damage					PASS	
Check TEE control, Tip Deflection					PASS	
Verify TEE Handle & Connector for Fluid					PASS	
Verify Control Knob Locks					PASS	
Review TEE Care & Handling with Customer					PASS	

Min	Max	UoM	Measurement	Result	Comments
-----	-----	-----	-------------	--------	----------

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## Debrief

Number of Soft Failures 2

Philips Representative Signature:

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